

# VIRACLEAN®

HOSPITAL GRADE DISINFECTANT  
AUST R 69000

## **ABOUT VIRACLEAN:**

VIRACLEAN® is a major development in advanced cleaning and disinfecting technology from Whiteley Medical. It is the result of years of intensive research and development. A considerable investment in both time and money has been made in the interests of improving the safety and efficacy of hospital grade disinfectants. VIRACLEAN® can be used for both cleaning and disinfecting, demonstrating beneficial versatility for healthcare and associated environs.

VIRACLEAN® has been proven to kill a broad range of bacteria and viruses including Enterococcus Faecalis (VRE), MRSA, Hepatitis B Virus, Herpes Simplex Virus and the Influenza Virus. The revolutionary hepadnacial test method was developed at the University of Sydney and sponsored with funding from Whiteley Medical.

## **MICROBIOLOGICALLY TESTED:**

VIRACLEAN® has been proven to kill:

- Staphylococcus aureus (MRSA or Golden Staph)
- Pseudomonas aeruginosa
- E coli
- Enterococcus Faecalis (VRE)
- Acinetobacter
- Acetobacter
- Salmonella choleraesuis
- Proteus vulgaris
- Influenza virus
- Herpes Simplex virus
- Hepatitis B Group virus
- Candida Albicans

*\*Independent test reports available on request.*

## **PRODUCT DESCRIPTION:**

VIRACLEAN® is a pink liquid with a lemon fragrance. It is easy to use and is conveniently packaged in squeeze or spray bottles, 5 litre containers, and 15 litre containers. This product is not intended for use on therapeutic devices. INSTRUMAX® PINK is available as an instrument disinfectant low level for use in disinfecting non critical medical devices.

Also available: V-Wipes Hospital Grade Disinfectant Wipes in Flatpacks and Canisters.

Viraclean® - June 2015 1

**PRODUCT FEATURES:**

**VIRACLEAN®:**

- is non flammable.
- is suitable for use on healthcare surfaces.
- does not contain chlorine or oxidising biocides.
- will not cause staining on fabrics or carpets.
- has a neutral pH.
- is friendly to metals (non corrosive).
- can be used neat for use as a hospital grade disinfectant.
- has a mild lemon fragrance.

Some of this data has been internationally published.<sup>1</sup> VIRACLEAN® passes the TGA option B (dirty conditions). This independent data is available on request.

**APPLICATION:**

The use of gloves is recommended whenever using disinfectants.

**TO CLEAN AND DISINFECT SURFACES:**

Use undiluted for disinfection.

To clean and remove organic soiling (blood, sweat, sputum etc) from surface, apply liberally around the soiled area and work towards the centre. Soak up using fresh, disposable cloth. Once clean, to disinfect the surface, apply undiluted and leave to dry for 10 minutes at room temperature.

If required, soak up using fresh, disposable towelling. Used towelling should then be disposed to a suitable receptacle.

**CHEMICAL SPECIFICATIONS:**

Surfactant Blend.....	Blended
Colour.....	Light Pink
pH.....	Neutral (7)
Specific Gravity.....	0.99
Corrosion Inhibitor.....	Not necessary
Foaming Properties.....	Slight
Occupational Safety.....	Use as spray or liquid
Dangerous Goods Rating.....	Not applicable
Shelf Life.....	18 months
Safety Equipment.....	Use of gloves and eye protection recommended

**VERSATILE:**

VIRACLEAN® is formulated to replace most cleaning products and hospital grade disinfectants. VIRACLEAN® can be used for cleaning and disinfecting.

Viraclean® - June 2015 2

**HELPS CONTROL INFECTION:**

VIRACLEAN® is the only solution that can be used for routine cleaning and can also be used for disinfecting where necessary. This makes for simplification of cleaning and disinfecting procedures.

**EDUCATIONAL VIDEO AVAILABLE:**

An Educational video is available on 'How to clean and disinfect blood spills using VIRACLEAN®'. To watch the educational video visit [www.whiteley.com.au/educational\\_videos](http://www.whiteley.com.au/educational_videos)

**PRODUCT USE:**

<i>Nursery:</i>	Cleaning cribs, cots, change tables.
<i>CSSD:</i>	Bench tops, stainless steel surfaces, sinks, trolleys, tap fittings.
<i>Toilet Shower:</i>	Sinks, tile surfaces, and all hard surfaces.
<i>Intensive Care Units:</i>	Cleaning and disinfection of all hard surfaces.
<i>Manual Cleaning of Theatre Equipment:</i>	VIRACLEAN® can be used to clean and disinfect equipment and cables and will not adversely affect surfaces.
<i>Operating Theatres:</i>	Cleaning and disinfecting of tables in theatre.
<i>Blood Spills:</i>	VIRACLEAN® is ideal for cleaning and disinfecting blood spills on all hard surfaces.
<i>Telephones:</i>	VIRACLEAN® is ideal for cleaning and disinfecting phones where large numbers of bacteria are found.
<i>Dental &amp; GP Surgeries:</i>	VIRACLEAN® can be used in Dental & GP Surgeries for cleaning and disinfecting surfaces including chairs, stools, tables, benches and light fittings.
<i>Laboratories &amp; Pathology Clinics:</i>	VIRACLEAN® is excellent for wiping down benches, disinfection of white rooms, cleaning and disinfection of fume cupboards, etc.

**SAFETY DATA:** Information is listed in the Material Safety Data Sheet.

**PACK SIZE:** 12 x 500mL Squeeze, 12 x 500mL Spray, 2 x 5 Litre, 15 Litre.

**ENVIRONMENTAL STATEMENT:**

Whiteley Corporation manufactures products that comply with Australian and International environmental policy. Whiteley Corporation aims to ensure that every product is 100% biodegradable, whilst achieving optimal performance and benefit to customers. All Whiteley Corporation products are supplied in 100% recyclable containers. Whiteley Corporation and its staff are fully committed to ensuring that our Australian environment and those international environs in which our products are used are enhanced by the use of technology, and that these environments are preserved for the enjoyment of future generations.