

Terms and Conditions of Labour TENS Hire

1. The Labour TENS machine will only be used by the person named on the hire form
2. The hirer will have discussed and sought approval from their medical practitioner before using the Labour TENS
3. The hirer will follow the instructions as provided and will use the machine under medical supervision
4. The success of Labour TENS cannot be guaranteed
5. **Return of the Labour TENS is the responsibility of the hirer.** A prepaid registered post return satchel is provided for your protection. Liability transfers from the hirer to Australia Post only when the return satchel is 'processed' by Australia Post. Should Labourtens.com.au not receive the returned TENS and Australia Post have no record of the item then the hirer shall be required to pay a replacement fee of up to \$200
6. Should the Labour TENS be lost or stolen you will be charged a replacement fee of up to \$200
7. Should the Labour TENS be returned damaged, contaminated with body fluids or with missing parts Labourtens.com.au reserve the right to charge additional fees for cleaning, repairs or replacements
8. Please advise the office as soon as possible if you require extra time. Additional fees of \$15 per week will apply if the Labour TENS is returned late.
9. In the event of a booking cancellation you will be refunded less an administration fee of \$10. **If the Labour TENS has already been dispatched or if the machine is unused NO REFUND can be made**
10. LabourTENS.com.au are not accountable for any injury or damage sustained as a result of misuse

Precautions and Contra-indications

1. Labour TENS contains electrical components and **MUST NOT** be used in **water** Remove TENS before entering bath or shower.
2. **SHOULD NOT** be used by persons with a **cardiac pacemaker**
3. Persons with **heart disease** or **epilepsy** should seek approval from their doctor before use
4. **SHOULD NOT** be applied internally or used on **broken** or **cracked skin**.
5. **SHOULD NOT** be placed on areas of the skin where normal sensation is absent
6. Pads **SHOULD NOT** be placed over the heart or chest, on the neck or head
7. Keep out of reach of children and pets
8. Electrical foetal monitoring may be affected by Labour TENS use. Please turn off TENS if requested by your midwife or doctor
9. **NOT** to be used before 35 weeks gestation without prior consultation from your doctor (other than short test)